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- 3) Do not use if you have any of the following conditions: diabetes, tuberculosis, benign or malignant tumors, phlebitis or thrombosis, hemorrhages, open or fresh wounds, ulcerated sores, varicose veins, poor circulation or bruised, discolored, burned, broken, swollen or inflamed skin.
- 4) This product should not be used by any individual suffering from a condition that limits or altogether inhibits the user’s capacity to feel or have sensation in any part of the body.
- 5) If you experience pain in a muscle or joint for a prolonged period of time, discontinue use and consult your doctor. Persistent pain could be a symptom of a more serious condition.
- 6) Use of this product should be pleasant and comfortable. Should pain or discomfort result, discontinue use and consult your physician.
- 7) Never use while sleeping or drowsy.
- 8) Ensure that your hands are dry when operating the controls or removing the plug.
- 9) Do not use while bathing or showering. Use only in a dry area.
- 10) Always unplug the unit from the electrical outlet immediately after using and before cleaning.
- 11) Keep out of reach of children.

### MAINTENANCE INSTRUCTIONS

- 1) Store massager in a secure place away from excessive heat or sharp objects. Do not handle the massager by any of the electrical cords and do not wrap the cord around the unit.
- 2) To clean the unit, make sure the massager is unplugged. Allow the unit to cool and then wipe down with a damp cloth. Do not use any strong chemical cleaners such as benzine.

### RECOMMENDED MASSAGING TIMES

Generally, massage should last 15 minutes on any one particular part of the body. Treatment may be applied once or twice daily as desired.

As metabolism varies among individuals, so does their response to massage. You will soon adjust the time to suit your particular needs and personal comfort. If you have any doubts regarding treatment, we recommend you consult your physician.

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## SAVE THESE INSTRUCTIONS

### OPERATING INSTRUCTIONS

To avoid risk of electric shock this appliance has been manufactured with a polarized plug (one blade is larger than the other). If the plug does not fit properly into your outlet do not force it or change in any way (seek a professional electrician to install the proper outlet). A 6-foot cord is provided for convenience.

Your deep tissue foot massager has been designed to bring you soothing comfort from your feet up.

Prior to use make sure your massager is placed on the floor or a level surface where added weight can be applied.

Plug unit into power source (120V A.C.).

Remove your shoes; your foot massager may be used with or with out socks depending on your preference.

If you desire heat with your massage, press the heat button first. Red infrared lights will indicate that the unit is heating up.

For soothing massage press the vibration button, one press for high, two presses is for low. Buttons are located at the base of the unit.

### CAUTIONS

- 1) Consult your doctor before using in case of illness or medical condition, including but not limited to pacemaker use, pregnancy, cancer, infection, fracture or persistent pain.
- 2) Do not use on any unexplained pain or swollen muscles, or following a serious injury, before consulting your physician.

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## IMPORTANT SAFETY INSTRUCTIONS

When using an electrical appliance, basic precautions should always be followed, including the following:

Read all instructions before using this massager.

## DANGER—

To reduce the risk of electric shock:

- 1) Always unplug this appliance from the electrical outlet immediately after using and before cleaning.
- 2) Do not reach for an appliance that has fallen into water. Unplug immediately.
- 3) Do not use while bathing or in a shower.
- 4) Do not place or store appliance where it can fall or be pulled into a tub or sink. Do not place in or drop into water or other liquid.

## WARNING—

To reduce the risk of burns, fire, electric shock, or injury to persons:

- 1) ⚠ NOT stand on the unit.
- 2) ⚠ NOT put sharp objects on or near mesh.
- 3) An appliance should never be left unattended when plugged in. Unplug from outlet when not in use, and before putting on or taking off parts.
- 4) Do not operate under blanket or pillow. Excessive heating can occur and cause fire, electric shock, or injury to persons.
- 5) This appliance should not be used by, on, or near children, or individuals with certain disabilities.
- 6) Use this appliance only for its intended use as described in this manual. Do not use attachments not recommended by Conair.
- 7) Never operate this appliance if it has a damaged cord or plug, if it is not working properly, or if it has been dropped or damaged, or dropped into water. Return the appliance to a service center for examination and repair.
- 8) Do not carry this appliance by supply cord or use cord as a handle.

CONAIR

BODY BENEFITS

deep tissue foot massager with microban®

INSTRUCTION BOOKLET FOR MODEL FM18

For your safety and continued enjoyment of this product, always read the instruction book carefully before using.

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Malla flexible para  
mayor comodidad

LIMITED ONE YEAR WARRANTY

Conair will repair or replace (at our option) your unit free of charge for 12 months from the date of purchase if the appliance is defective in workmanship or materials. To obtain service under this warranty, return the defective product to the service center listed below, together with your purchase receipt and \$9.95 for postage and handling. California residents need only provide proof of purchase and should call 1-800-326-6247 for shipping instructions. In the absence of any implied warranties, obligations or liabilities, including but not limited to the implied warranty of merchantability and fitness for a particular purpose, shall be limited in duration to the 12-month duration of this warranty. Some states do not allow limitations on how long an implied warranty lasts, so the above limitations may not apply to you.

IN NO EVENT SHALL CONAIR BE LIABLE FOR ANY SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES FOR BREACH OF THIS OR ANY OTHER WARRANTY, EXPRESS OR IMPLIED, WHATSOEVER. Some states do not allow the exclusion or limitation of special, incidental or consequential damages, so the above limitation may not apply to you.

This warranty gives you specific legal rights, and you may also have other rights which vary from state to state.

us

Service Center

Conair Corporation  
Service Department  
475 North Glen Harbor Blvd.  
Glendale, AZ 85307

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us

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Conair Corporation  
Service Department  
475 North Glen Harbor Blvd.  
Glendale, AZ 85307

©2008 Conair Corporation

Visítenos en la Web en [www.conair.com](http://www.conair.com)

FM18  
08NP143

What is Microban®

Microban® antimicrobial protection is built-in during manufacturing to provide continuous antimicrobial protection. Your deep tissue foot massager has been infused at the footpads with Microban. Microban resists bacterial odors, mold and mildew.

HEAT

Therapeutic infrared heat soothes the body from your feet up

Heat switch

Vibration switch

Flexible mesh for added comfort

MESH

Massage

Soft, comfortable gel pads

Model FM18

DEEP TISSUE FOOT MASSAGER WITH MICROBAN®

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